





K9 HELM™




Place template on snout,
about 6 centimeters
in front of eyes. Try each size
and decide which one fits best.



Template #

S1


Snout Width



Template #

S2

Snout Width




IMPORTANT

Template #

S3

Snout Width



This measurement must be verified with a ruler after printing!
PDF must be printed at 100% with NO PAGE SCALING to fit paper size.

Cut along dotted line and use upper portion as shown in illustration above. For best results, after cutting to shape, trace onto cardboard and cut out. Be sure to write the template numbers on the cardboard cutout.